



# Calibration


Nicholas W. M. Ritchie  
Physicist, Microanalysis Research Group  
National Institute of Standards and Technology  
Gaithersburg, MD 20899-8371  
*[nicholas.ritchie@nist.gov](mailto:nicholas.ritchie@nist.gov)*



An x-ray spectrum is actually an accumulation of thousands to millions of individual x-ray event presented as a histogram. The x-axis of the histogram is a series of bins in which a number of x-ray events are accumulated. The width of each bin is typically nominally either 5 eV or 10 eV.



X-rays go into bins dependent on their energy...



# Calibrated parameters

- The width of each bin (“channel width”)
  - The width is typically expressed in eV/bin
- The offset of the 0<sup>th</sup> bin (“zero offset”)
  - The offset is typically expressed in eV.
  - Many detectors have a negative offset of a few hundred eV. Zero offset is also common.
- The resolution of the detector at Mn K $\alpha$  (“FWHM”).
  - The resolution is generally reported as the full-width-at-half-maximum (background corrected) for the Mn K $\alpha$  peak.

These parameters can be extracted from a suitable measured spectrum

# Calibrating what?

## **Calibrating the detector**

EDS detector vendors usually provide a mechanism for adjusting the hardware to produce a bin width close to the nominal value.

On modern detectors the calibration is often performed automatically in the vendor software.

Although DTSA-II's detector model can compensate for a poorly calibrated detector, in practice, the detector should always be calibrated using the vendor's software.

## **Calibrating the DTSA-II's detector model**

To interpret and simulate spectra, DTSA-II needs to know the real performance specifications of the detector. The calibration tool extracts this information from a measured spectrum.

File Process Tools Report Help

Quantification alien  
Simulation alien  
Calibration alien  
Optimization alien  
Quality control alien  
Edit spectrum properties  
Assign material  
Edit standards database

Counts

5 10  
15 201000  
0  
2  
4

DTSA-II provides a tool for  
extracting the detector  
calibration from a measured  
spectrum.

BaF2 15.0 keV 1.000 nA 0

12

14

Spectrum Report Command

Default Detector

TESCAN MIRA3

SDD (Medium, 4096)

Spectrum List

BaF2 15.0 keV 1.000 nA 0

Spectrum Properties

Name

Acquisition time

5/28/14 8:20 AM

Aluminum layer thickness

0 nm

Aluminum window thickness

30 nm

Azimuthal angle

0°

Beam energy

15.0 keV

Calibration GUID

caf24d30-9fad-62b9-2f76-dac1a64cbf74

Dead layer

0 µm

Detector

SDD (Medium, 4096) - FWHM[Mn Kα] = 130.8...

Detector GUID

b8b98369-7eb7-32eb-d2c8-c9d91a7e78f8

Detector area

30 mm²

Detector orientation

[-0.819,-0.000,0.574]

Detector position

[27.851,0.000,-2.502]

Detector thickness

0.45 mm

Detector type

Silicon Drift Detector

Detector window

Moxtek AP 3.3 (manufacturer's table)

Display name

BaF2 15.0 keV 1.000 nA 0

Duane-Hunt

15.646 keV

Elevation

35°

Energy offset

-0.4 eV

KLM Lines

Element

F

☒ F☒ K-family☐ Edges☒ Z-order☐ E-order

Fluorine

Clear

Clear All

Standard Composition

Element	Mass Fraction	Atomic Fraction
Fluorine	0.2167	0.6667
Barium	0.7833	0.3333

1 spectra selected

The "calibration alien" will step you through the process of calibrating a detector.

## Calibrate an EDS detector

First page



## Select a detector

Next: Calibration method

## Instrument and Detector

Instrument: TESCAN MIRA3  
Detector: SDD (Medium, 4096) - FWHM[Mn Ka]=130.8 eV - initial

Message: Select a detector to calibrate

First, specify which detector on which instrument you wish to calibrate.

Spectrum Report Command

Default Detector


TESCAN MIRA3

SDD (Medium, 4096)

Spectrum List

BaF2 15.0 keV 1.000 nA 0

1 spectra selected





NIST DTSA-II

File Process Tools Report Help

Calibrate an EDS detector

Previous: Calibration

Next: A

**Measured spectrum**

Specify a spectrum

Spectrum Pure copper[Wed Apr 23 080417 2014][0] ...

Specify the material

Material Cu standard

Live time 59.24 sec. Probe current 0.715 nA

Fit type Linear (default)

Specify an effective date

Effective date Apr 23, 2014

Message: The reference beam energy is 20.0 keV.

More..

Back Next Finish

1 spectra selected

5000  
4000  
3000  
bunts

893 894

BaF2 15.0 keV 1.000 nA 0

12 14

SDD (Medium, 4096)

Spectrum List

BaF2 15.0 keV 1.000 nA 0

Leave the "Fit type" as "Linear" unless you have a good reason to do otherwise.

The composition of the material from which the spectrum was collected. You can select one of the default materials or enter an arbitrary material.

A spectrum file collected from a material of known composition. Transition metals are ideal calibration standards.

This information should be extracted from the spectrum file or provided by you.

The date on which the calibration will become effective. Read from the spectrum file or set manually.

BaF2 15.0 keV 1.000 nA 0  
Pure copper[Wed Apr 23 080417 2014][0]  
BestFit[Pure copper[Wed Apr 23 080417 2014][0]]  
Characteristic[Pure copper[Wed Apr 23 080417 2014][0]]

Calibrate an EDS detector

Previous: Fit Results

## Fit results

Next: Spec

This detector has nominally a 0 eV offset and a 10 eV/bin channel width. You should expect the measured channel width to be within a few thousands of the nominal value.

## Energy calibration

Fano factor 0.1200

Noise 6.32±0.05 eV

FWHM 131.25±0.17 eV

## Resolution calibration

Zero offset 1.61±0.29 eV

Channel width 9.9996±0.0005 eV/channel

☒ Add to database☐ Output elemental fits

More...

Back

Next

Finish

Cancel

Silicon Drift Detector  
Moxtek AP 3.3 (manufacturer's table)  
35°  
-0.4 eV  
10 eV/channel  
0 nm  
TESCAN MIRA3

The algorithm will fit the Bremsstrahlung and characteristic peaks and report summary information here.

The FWHM is reported as though the material was Mn. The algorithm actually fits the "noise" and calculates the FWHM from the noise and the Fano factor.

NIST DTSA-II

File Process Tools Report Help

50 000 Same scale

Calibrate an EDS detector

Results

Energy calibration

Fano factor 0.1200

Noise 6.32±0.05 eV

FWHM 131.25±0.17 eV

Resolution calibration

Zero offset 1.61±0.29 eV

Channel width 9.9996±0.0005 eV/channel

☒ Add to database

☐ Output elemental fits

More...

Back Next Finish Cancel

5 spectra selected

Default Detector

TESCAN MIRA3

SDD (Medium, 4096)

Spectrum List

BaF2 15.0 keV 1.000 nA 0

Pure copper [Wed Apr 23 08:04:17 2014]

Bremsstrahlung [Pure copper [Wed Apr 23 08:04:17 2014]]

Characteristic [Pure copper [Wed Apr 23 08:04:17 2014]]

BestFit [Pure copper [Wed Apr 23 08:04:17 2014]]

Message:

More...


Mass Fraction Atomic Fraction

5 spectra selected

If the summary results seem reasonable, you probably want to add the calibration to the database.

By default, the raw spectrum, the modeled best fit, the characteristic spectrum and the Bremsstrahlung spectrum are added to the Spectrum List. If the material has more than one element and you would like to see the contribution of each element separately, check this box.

The "Add to database" check box is disabled if a non-linear bin to energy model is selected.




# What is the “Effective date?”

- Each data spectrum is associated with the most recent prior calibration as it is read in.
- The assumption is that the calibration spectrum is collected before or at least on the same day as the data spectra. You can fake this by manually changing the effective date.
- A data spectrum will be associated with the most recent calibration.
- Calibrations are pre-dated to take effect at 12:01 AM on their “effective date.” Each subsequent calibration collected on a day is assumed to take effect 1 hour later than the previous one.

# Effective date

- △ Calibration spectrum
- ⊙ Data spectrum



The first calibration on a day becomes effective at 12:01 AM  
The second becomes effective at 1:01 AM, third at 2:01 AM etc.

# How is the calibration used?

- Quantification
  - To determine the position and width of the region to fit for each characteristic x-ray transition.
- Simulation
  - To determine the bins into which the simulated Bremsstrahlung and characteristic x-rays will be placed.
  - To determine the simulated width of the characteristic x-ray peaks.
- KLM lines
  - Ensures the peaks line up with the KLM lines

# Advanced technique

The “calibration alien” “calibrate using standard spectrum” tool can actually fit spectra from arbitrary materials.

Use 1: This can be useful when a standard requires a reference but a good reference can not be found. You can use the “output elemental fits” to export “characteristic-only spectra” for any element which can be used as either a reference or a standard.

Use 2: The elemental fits can also be used to subtract off undesired characteristic peaks.