


# An Introduction to “Standard Bundles” in NIST DTSA-II

Nicholas W. M. Ritchie  
Microanalysis Group  
National Institute of Standards and Technology  
Gaithersburg, MD 20899  
*nicholas.ritchie@nist.gov*

Edited: 4-May-2016

# Using a complex standard

DTSA-II has added support for “standard bundles”, collections of spectra designed to simplify the process of quantifying unknown spectra using “complex” or “similar standards”



Eh? What did you say?

# Standards


- Spectrum quantification is based on comparing the *unknown spectra* with spectra from materials of *known composition*.
- To a good approximation, the composition in the unknown equals<sup>1</sup> the composition in the standard times the ratio of the intensity<sup>2</sup> in the unknown relative to the standard.

$$C_u \approx C_s \frac{I_u}{I_s}$$

<sup>1</sup> Ignoring “matrix corrections”

<sup>2</sup> Where *intensity* is understood to mean the intensity in a characteristic x-ray peak associated an element.


# Example – Simple standards



Element	Z	Mass Frac
Oxygen	8	0.34
Magnesium	12	0.0302
Silicon	14	0.187
Titanium	22	0.06
Zinc	30	0.0402
Zirconium	40	0.074
Barium	56	0.2687
Total		1.0001

## Standards:

Si	SiO <sub>2</sub>
Mg	MgO
Ti	Ti metal
Zn	Zn metal
Zr	Zr metal
Ba	Ba(Si <sub>2</sub> O <sub>5</sub> )
O	calculated



While SiO<sub>2</sub>, MgO and Ba(Si<sub>2</sub>O<sub>5</sub>) contain multiple elements, the elements within each standard do not interfere with each other so the peak shapes are equivalent to the pure element peak shapes.



## Select a quantification mode

First page

Next: *Specify the instrument*

Select the mode which best describes the operation you wish to perform. The mode you select will determine what information you will be asked to provide and what information will be computed.

- ☒ Determine the composition of an 'unknown' spectrum by MLLSQ fitting to standards
- ☐ Determine the composition from k-ratios
- ☐ Quantify a STEM spectrum using MLLSQ fitting and  $\zeta$ -factors

Message: Select an analysis mode.

More...

Back

Next

Finish

Cancel

Select quantify a bulk material using standards fit using multiple linear least squares fitting.

Quantification Alien

Quantification Alien

Previous: *Select a quantification mode*

Next: *Specify standard spectra*

## Specify the instrument

Instrument

Acquired on the

Detector

using the


with calibration

Setting

at a beam energy of  keV.

Message:

Verify that the default instrument and detector information read from the spectrum/spectra is correct.




Specify standard bundles using this button.

Alternatively, specify raw spectra using this button.


The standards you have specified show up here along with information about the standard.

You may be asked to provide compositional information if this information isn't present in the spectrum.




You may either choose to provide a standard for O  
or you may choose to have the quantity of O  
computed using an assumed stoichiometry (which  
you can edit if the default assumptions aren't right.)






Because we chose simple standards, there is no need to specify separate reference spectra. (The standardized elements don't interfere with other elements in the standard.)





You can allow DTSA-II to select which line to use to quantify, or you can force it to use the line that you specify.



This panel gives you an opportunity to ensure that all the necessary data is available for the unknowns.

Update spectrum properties here.






Name	O	Mg	Si	Ti	Zn	Zr	Ba
K240 O	32.0	2.7	17.2	5.9	4.0	7.1	27.4
K240 I	32.5	2.7	17.5	6.1	4.0	7.3	27.6
Average	32.3	2.7	17.3	6.0	4.0	7.2	27.5
Nominal	34.0	3.0	18.7	6.0	4.0	7.4	26.9
Difference	-5%	-11%	-7%	0%	0%	-3%	2%

$$C_u \approx C_s \frac{I_u}{I_s}$$

- The validity of this approximation is enhanced when the composition of the unknown and the standard are approximately equal.
- This condition is called a “similar standard”
- Similar standards don't need to contain all the elements in the unknown.
- Similar standards can contain elements not in the unknown.

# Example – Similar/Complex Standards



## Remotivate:

Selecting a standard with composition similar to our unknown reduces the need for matrix correction and is more likely to produce accurate results.

Like our unknown glass, Benitoite contains Ba and Ti.


	Benitoite	“Unknown”
Barium	33.22 %	26.9 %
Titanium	11.58 %	6.0 %

Not in exactly the same proportions but similar.

But there is a nasty overlap between Ba and Ti! Can we really use Benitoite as a standard?

# References

- Can we use Benitoite as a standard?
  - Yes, if we tell the program what pure Ba and pure Ti look like!



We provide “references” (Ti metal and  $\text{Ba}(\text{Si}_3\text{O}_9)$ ) that tell the program what unobstructed Ti K and unobstructed Ba L lines look like.





First page

## Select a quantification mode

*Next: Specify the instrument*

Select the mode which best describes the operation you wish to perform. The mode you select will determine what information you will be asked to provide and what information will be computed.

- ☒ Determine the composition of an 'unknown' spectrum by MLLSQ fitting to standards
- ☐ Determine the composition from k-ratios
- ☐ Quantify a STEM spectrum using MLLSQ fitting and  $\zeta$ -factors

Message: Select an analysis mode.

[More...](#)[Back](#)[Next](#)[Finish](#)[Cancel](#)

Same as before.



## Specify the instrument

Previous: *Select a quantification mode*

Next: *Specify standard spectra*

Instrument

Acquired on the MIRA-3

Detector

using the SDD (Medium, 4096)

with calibration FWHM[Mn K $\alpha$ ] = 130.6 eV - 2016-04-22 12:00

Setting

at a beam energy of 15.0 keV.

Message:

More...


Back

Next


Finish

Cancel


Same as before.




Same as before.




This time we specify Benitoite as a standard for Ba, Ti, Si and O. (Previously we'd compute O.) Mg is now pure Mg and Zn and Zr are also pure elements.




Now we're measuring O because we can measure O accurately with a similar standard.



This time we need to specify spectra as “peak shape” references to help us to extract the intensity information from both the standard and the unknown for Ba and Ti.





Same  
as before.






# Results

## Simple standards

Name	O	Mg	Si	Ti	Zn	Zr	Ba
K240 O	32.0	2.7	17.2	5.9	4.0	7.1	27.4
K240 I	32.5	2.7	17.5	6.1	4.0	7.3	27.6
Average	32.3	2.7	17.3	6.0	4.0	7.2	27.5
Nominal	34.0	3.0	18.7	6.0	4.0	7.4	26.9
Difference	-5%	-11%	-7%	0%	0%	-3%	2%

## Similar standards

Name	O	Mg	Si	Ti	Zn	Zr	Ba
K240 O	33.6	3.0	18.7	6.0	4.1	7.2	27.1
K240 I	34.6	3.0	19.0	6.2	4.0	7.3	27.3
Average	34.1	3.0	18.9	6.1	4.0	7.2	27.2
Nominal	34.0	3.0	18.7	6.0	4.0	7.4	26.9
Difference	0%	0%	1%	1%	0%	-2%	1%



WEE HEE!  
WOO HOO!

As is often the case, similar standards produce more accurate results than simple standards.

# Making similar standards easy

## The problem:


Similar standards often require reference spectra and specifying reference spectra every time is tedious!

## The solution:

Bundle standard spectra with the reference spectra required to use them.

## One or more spectra

A spectrum bundle can be created from one or more spectra collected under identical conditions from the same material. Highlight the spectra in the spectrum list before selecting "Make 'Standard Bundle'".



# A word about collecting standards


- Standards have to be high quality spectra collected under well controlled conditions from a homogeneous material.
- Many potential pitfalls can be detected by collecting multiple spectra from different points on the same material, comparing the spectra and later combining them into a single spectrum.
- It is better to collect 5 spectra of 60 seconds each, compare and combine them to form one 300 second standard than to collect one 300 second spectrum.
  - If all five spectra are essentially identical, it raises your confidence in the reliability of the spectra (no contamination, geometric effects, pits, ridges, current drift,...)
  - Even if you end up discarding one or two spectra, you feel confident the other three or four are of high quality.

# Describing the standard

Specify the composition of the standard.

Specify the elements for which this material can be a standard.

Also verifies that all the necessary instrument data is present.




# Describing the standard

Specify the composition of the standard.

Specify the elements for which this material can be a standard.

Also verifies that all the necessary instrument data is present.




# Specifying references

Specifies the references necessary to use this standard.

Specifies the reference spectra for each “region-of-interest.”


Allows you to also strip spurious C or O if present in the standard spectrum.



# Specifying references

Auto-create references by fitting an analytical model to the standard.

Read the references from spectrum files on disk.






# Specifying references

Auto-create references by fitting an analytical model to the standard.

Read the references from spectrum files on disk.




# Saving the standard bundles

When you are done, save the standard bundles to disk.

You will save one bundle of each element for which the spectrum is a standard.

By default, the bundles are named by element and material.



# To reiterate

- Standard bundles are designed to simplify using the “Quantification Alien”
- Standard bundles contain:
  - A standard spectrum
  - Zero or more references spectra necessary to use the standard
- References may be missing from the bundle
  - In this case, you will be asked to provide the missing reference when you use the standard bundle in the “Quantification Alien”
- Standard bundles are so convenient that it is generally easiest (but not strictly necessary) if all standard spectra (simple or complex) are transformed into standard bundles.