

# **STANDARDISATION and TRACEABILITY in LABORATORY MEDICINE**



**Mathias M. Müller**  
**Institute of Laboratory Diagnostics**  
**Vienna - Austria**

# **LABORATORY MEDICINE**

**... an integrated discipline in health care:  
risk assessment**

**diagnosis of health and disease,  
follow-up and monitoring of patients.**

**... using physical, chemical, biochemical,  
immunological, molecular biological techniques  
for measurements of  
body fluids, tissues, and cells**

# **ANALYTES - MEASURANDS**

## **CATEGORY I**

**Traceable to the SI**  
**Reference Systems**

**creatinine**

**cholesterol**

**total/fractions**

**enzymes**

**electrolytes**

**total/activity/free**

**glucose**

**steroids & thyroxine**

**free/bound to proteins**

## **CATEGORY II**

**Not-traceable to the SI**  
**Conventional RMs**

**(international agreement)**

**coagulation factors**

**glycoproteins**

**isoforms/glycoforms**

**proteins**

**peptide-bond (biuret reaction)**

**epitopes (antibodies)**

# BIOLOGICAL VARIATIONS

## **Matrix - Fluids**

Serum - Plasma

Urine

Liquor

Ascites

Tissue

Cells

## **Factors**

Age

sex

time of day (circadian rhythm)

posture (supine/upright)

serum/plasma

fasting/non-fasting


weight/overweight

seasonal influence

# **Intra-individual variations**


	<b>Analytical CV</b>	<b>Daily %</b>	<b>Weekly %</b>	<b>Monthly %</b>
<b>Na</b>	<b>0.6</b>	<b>1.4</b>	<b>0.8</b>	<b>1.3</b>
<b>K</b>	<b>1.0</b>	<b>7.8</b>	<b>6.7</b>	<b>7.3</b>
<b>Glucose</b>	<b>1.5</b>	<b>25.8</b>	<b>16.8</b>	<b>20.8</b>
<b>Creatinine</b>	<b>1.6</b>	<b>6.8</b>	<b>6.9</b>	<b>13.6</b>
<b>Uric acid</b>	<b>1.0</b>	<b>9.8</b>	<b>12.4</b>	<b>14.3</b>
<b>ALT (enzyme)</b>	<b>0.9</b>	<b>10.3</b>	<b>32.2</b>	<b>47.5</b>

# ANALYTICAL BIAS



IMEP- 17: Trace and minor constituents in human serum

Certified value :  $140.36 \pm 0.95 \text{ mmol} \cdot \text{L}^{-1}$  [ $U = k \cdot u_c$  ( $k=2$ )]




Results from all participants (992 laboratories)

# Amylase Comparison

IMEP- 17: Trace and minor constituents in human serum

Certified value :  $56.8 \pm 2.6 \text{ U} \cdot \text{L}^{-1}$  [ $U = k \cdot u_c$  ( $k=2$ )]



All reported results (863) arranged in method groups:

IFCC comparable methods; Different methods, Scandinavian level;

Different methods, Roche level; Different methods, Original level; Vitros 250-950, calculated to

IFCC; Vitros 250-950, Scandinavian level; Vitros 250-950, calculated to Roche level;


Vitros 250-950, original level and Other/No info



# Amylase Comparison

IMEP- 17: Trace and minor constituents in human serum


Certified value :  $56.8 \pm 2.6 \text{ U} \cdot \text{L}^{-1}$  [ $U=k \cdot u_c$  ( $k=2$ )]



Results from participants from EU Countries: Austria, Belgium, Denmark, Finland, Germany.

IMEP- 17: Trace and minor constituents in human serum


Certified value :  $56.8 \pm 2.6 \text{ U} \cdot \text{L}^{-1}$  [ $U=k \cdot u_c$  ( $k=2$ )]



Results from participants from North America; Canada and USA

IMEP- 17: Trace and minor constituents in human serum

Certified value :  $56.8 \pm 2.6 \text{ U} \cdot \text{L}^{-1}$  [ $U=k \cdot u_c$  ( $k=2$ )]



Results from participants from Nordic Countries;  
Denmark, Finland, Iceland, Norway and Sweden

- Difference in field methods
- Patient results not comparable




**Need for Harmonisation**

# Analytical Bias - Therapeutic Consequence

**HbA1c**


**Insulin-dependent  
patients**

**P. Hyltoft  
Petersen et al  
1997**




# Cholesterol Measurements

## Improved Cholesterol Measurement Accuracy Saves Health Care Dollars



Improvement in precision since 1968 has been estimated to save \$100M/yr in treatment costs

# Bias in Cholesterol Measurement Effects Medical Decision-Making




**Cholesterol Frequency  
Distribution of >20,000  
Mayo Clinic Patients**  
(with +1%, +3% and +10% limits  
around 240 mg/dL criteria point)

<u>If measurement bias were:</u>	<u>Positives (&gt;240 mg/dL) per 1000</u>	<u>Predicted Change in "Positives/1000"</u>
-10% bias	120	-129
-3% bias	203	
-1% bias	234	
0% bias	249	
+1% bias	263	+197
+3% bias	300	
+10% bias	446	

## AQAS - Method Target Values in 2 Control Samples

Analyte	LIA	EIA	MEIA	RIA
AFP (ng/ml)	39 79	32 93	35 99	26 72
CEA (ng/ml)	4 28	5 31	4 25	4 28
CA 19-9 (U/ml)	12 63	21 84	26 87	15 52
PSA (ng/ml)	12 120	7 70	7 80	

# Performance of 2 Field Methods



**Diagnostic procedures**

**Costs - €**

**CA 19-9 tests**

**246**

**X-Ray**

**157**

**Ultrasound**

**1062**

**Computer tomography**

**1062**

**NMR**

**1194**

**Consequences for the Patient**

**2,871**

# **EQUAS Results**

## **Clinical Guidelines for Decisions**



### **NEED FOR INTERNATIONAL STANDARDISATION**



- **Characterisation of Analyte**
- **Clinical Needs**
- **Reference Procedure**
- **Reference Material**
- **Reference Laboratories**

# STANDARDISATION

A technical process to reach conformity of  
measurement procedures by applying  
highest scientific standards

**REFERENCE SYSTEM**

**REFERENCE METHODS**  
**REFERENCE MATERIALS**  
**REFERENCE LABORATORIES**



- **ISO/EN 15195**

Requirements for **reference measurement laboratories** in laboratory medicine

- **EN 12286**

Measurements of quantities in samples of biological origins – Presentation of **reference measurement procedures**

- **EN 12287**

Description of **reference materials**



- **FULL NATIONAL MEMBER SOCIETIES: 79**
- **AFFILIATE MEMBER SOCIETIES: 4**
- **CORPORATE MEMBERS: 38**
  
- **REGIONAL ORAGNISATIONS  
AFFILIATED WITH IFCC**
  - ▶ **Arabic Federation of Clinical Biology**
  - ▶ **Asian Pacific Federation of Clinical Biochemistry**
  - ▶ **Colabiocli - Latin American Federation**
  - ▶ **FESCC - Federation of European Societies**

*ifcc*

# SCIENTIFIC DIVISION REFERENCE METHODS

- **ENZYMES**


- ALAT, ASAT, Amylase, AP, CK, gGT, LDH
- Lipase (in preparation)

- **BLOOD GASES - ELECTROLYTES**

- Tonometry
- $p\text{CO}_2$
- Na, K,  $\text{Ca}_{\text{ionised}}$

- **PROTEINS**

- Apo A1
- Hb
- HbA1c



# AQAS Interlaboratory Comparison

## Source: Austrian Clinical Chemistry Surveys

<b>ANALYTE</b>	<b>1970</b>	<b>1982</b>	<b>1992</b>	<b>2001</b>	<b>2003</b>
Participants	36	269	603	1358	1812
<b>AP</b>	33.9	11.5	5.6	5.9	5.1
<b>ALAT</b>	67.8	17.2	5.0	5.9	5.7
<b>ASAT</b>	36.3	20.8	5.3	6.0	5.8
<b>LDH</b>	32.0	10.1	5.9	3.7	4.0
CV %					




# SCIENTIFIC DIVISION

## REFERENCE MATERIALS

Apo A1	WHO: SP1	The logo of the World Health Organization (WHO), featuring a blue globe with a yellow caduceus in the center, surrounded by a laurel wreath.
Apo B	WHO: SP 3	
Albumin	WHO: 74/1	
Plasma Proteins	IRMM: CRM 470	The logo for the International Reference Materials Institute (IRMM), featuring the lowercase letters 'irm' in a stylized, italicized font.
PSA, free, complexed	WHO: 96/668, 96/670	
ALAT	IRMM: 454	
Amylase	IRMM: 456	A photograph of laboratory glassware, including a test tube and a vial, with a pipette tip visible, set against a dark background.
CK-MB	IRMM: 455	
gGT	IRMM: 452	
LDH-1	IRMM: 453	
Cortisol	IRMM: 451	
HCG primary standards	WHO: 99/642, 650, 688, 692, 708, 720,	
Lp(a)	WHO: 03/	
ASAT, HbA1c, myoglobin	IRMM: in preparation	


# CRM 470's EFFECT ON CALIBRATION



***ifcc SD***


## WG-CALIBRATORS IN CLINICAL ENZYMOMOLOGY (WG-CCE)

### Effect of Uniform Calibration on Various Methods



# TRACEABILITY

Number of Measurements



Property of the result related to national or international standards through an unbroken chain of comparisons all having stated uncertainties.

**A: traceable to SI**

**B: non-traceable to SI**

- Int'l Reference measurement procedure and int'l calibrator
- Int'l Reference measurement procedure but no int'l calibrator
- Int'l calibrator but no int'l reference measurement procedure
- Manufacturer's measurement procedure but neither int'l reference measurement procedure nor int'l calibrator



# IVD-Directive 98/79

The **traceability** of values assigned to calibrators and or control materials must be assured through **reference measurement procedures and reference materials** of a higher order

## ISO Standards

*In vitro* diagnostic medical devices - Measurements of quantities in biological samples

- **ISO 17511**– **Metrological traceability** assigned to calibrators and control materials.
- **ISO 18153** – **Metrological traceability** of values for catalytic concentration of enzymes assigned to calibrators and control materials.

# **JOINT COMMITTEE on TRACEABILITY in LABORATORY MEDICINE**

**JCTL**  
a global initiative,  
formed  
in Paris, June 12, 2002

[http://www.bipm.org/enus/2\\_Committees/JCTL.html](http://www.bipm.org/enus/2_Committees/JCTL.html)

# JCTLM

## Joint Committee for Traceability in Laboratory Medicine




Establishment of  
a global

**REFERENCE SYSTEM,**

**TRACEABILITY**



needs collaboration and mutual  
recognition between



Professionals, Metrology  
Institutes, Regulators, and IVD-  
Industry

**A JOINT VENTURE OF PROFESSIONALS**

# Focus on Standardisation and Traceability

- ◆ Excellence in Analytical Performance based on modern concepts of metrology and science
- ◆ Needs for Patients
- ◆ Impact on Clinical Decisions

...will add **QUALITY** and **VALUE** to  
**CLINICAL CHEMISTRY**  
and  
**LABORATORY MEDICINE**

