



Probabilistic Software Workshop  
September 29, 2014

Probabilistic Genotyping – An Overview  
Michael D. Coble, Ph.D.

---

---

---

---

---

---

---

---

### Product Disclaimer

- I will mention software programs and STR kit names and information, but I am in no way attempting to endorse any specific products.
- **NIST Disclaimer:** Certain commercial equipment, instruments, software programs, and materials are identified in order to specify experimental procedures as completely as possible. In no case does such identification imply a recommendation or it imply that any of the materials, instruments or equipment identified are necessarily the best available for the purpose.
- **Points of view are mine** and do not necessarily represent the official position of the National Institute of Standards and Technology or the U.S. Department of Justice. **Our group receives or has received funding from the FBI Laboratory and the National Institute of Justice.**

---

---

---

---

---

---

---

---

### Two Parts to Mixture Interpretation

- Determination of alleles present in the evidence and **deconvolution of mixture components** where possible
  - Many times through comparison to victim and suspect profiles
- **Providing some kind of statistical answer** regarding the weight of the evidence
  - There are multiple approaches and philosophies

---

---

---

---

---

---

---

---

### Statistical Approaches with Mixtures

See Ladd et al. (2001) Croat Med J. 42:244-246

<p><b>“Exclusionary” Approach</b></p> <p><b>Random Man Not Excluded (RMNE)</b></p> <p><i>Combined Prob. of Inclusion (CPI)</i></p> <p><i>Combined Prob. of Exclusion (CPE)</i></p>	<p><b>“Inferred Genotype” Approach</b></p> <p><b>Random Match Probability (RMP)</b></p> <p><b>(mRMP)</b></p> <p><b>Likelihood Ratio (LR)</b></p>
--	--

---

---

---


---

---

---

---

---




---

---

---

---

---

---

---

---



---

---

---

---

---

---

---

---



### The "2p" Rule

- "This rule arose during the VNTR era. At that time many smaller alleles "ran off the end of the gel" and were not visualised."

- Buckleton and Triggs (2006)

Is the 2p rule always conservative?"

---

---

---

---


---

---

---

---

### The "2p" Rule



Stain = AA

Suspect = AA

ST

LR = 100

$f(a) = 0.10$   $1/p^2 = 100$   $1/2p = 5$

---

---

---

---


---

---

---

---

### The "2p" Rule



Stain = AA

Suspect = AB

ST

Exclusion

$f(a) = 0.10$   $1/2p = 5$

---

---

---

---

---

---

---

---

# Likelihood Ratio

---

---

---

---

---

---

---

---

## Likelihood Ratios in Forensic DNA Work

- We evaluate the evidence ( $E$ ) relative to alternative pairs of hypotheses
- Usually these hypotheses are formulated as follows:
  - The probability of the evidence if the crime stain originated with the suspect or  $\Pr(E|S)$
  - The probability of the evidence if the crime stain originated from an unknown, unrelated individual or  $\Pr(E|U)$

$$LR = \frac{\Pr(E|S)}{\Pr(E|U)}$$

← The numerator  
← The denominator

Slide information from Peter Gill

---

---

---

---

---

---

---

---

## Likelihood Ratio (LR)

- Provides ability to express and evaluate both the prosecution hypothesis,  $H_p$  (the suspect is the perpetrator) and the defense hypothesis,  $H_d$  (an unknown individual with a matching profile is the perpetrator)

$$LR = \frac{H_p}{H_d}$$

- **The numerator,  $H_p$ , is usually 1** – since in theory the prosecution would only prosecute the suspect if they are 100% certain he/she is the perpetrator
- The denominator,  $H_d$ , is typically the profile frequency in a particular population (based on individual allele frequencies and assuming HWE) – i.e., **the random match probability**

Slide information from Peter Gill

---

---

---


---

---

---

---

---




---

---

---

---

---

---

---

---

Forensic Science International: Genetics 6 (2012) 679–688

DNA commission of the International Society of Forensic Genetics:  
 Recommendations on the evaluation of STR typing results that may  
 include drop-out and/or drop-in using probabilistic methods

P. Gill<sup>a,b\*</sup>, L. Gusmão<sup>c</sup>, H. Haned<sup>d</sup>, W.R. Mayr<sup>e</sup>, N. Morling<sup>f</sup>, W. Parson<sup>g</sup>, L. Prieto<sup>h</sup>,  
 M. Prinz<sup>i</sup>, H. Schneider<sup>j</sup>, P.M. Schneider<sup>k</sup>, B.S. Weir<sup>l</sup>

---

---

---

---

---

---

---

---

Summary of recommendations of the  
 ISFG DNA commission

- (1) Probabilistic methods following the 'basic model' described here can be used to evaluate the evidential weight of DNA results considering drop-out and/or drop-in.
- (2) Estimates of drop-out and drop-in probabilities should be based on validation studies that are representative of the method used.

---

---

---

---

---

---

---

---

Summary of recommendations of the  
ISFG DNA commission

- (3) The weight of the evidence should be expressed following likelihood ratio principles.
- (4) The use of appropriate software is highly recommended to avoid hand-calculation errors.

---

---

---

---

---

---

---

---

Probabilistic Approaches

- “Semi-Continuous” or “Fully Continuous”
- Semi-Continuous – information is determined from the alleles present – peak heights are not considered.
- Fully Continuous – incorporation of biological parameters (PHR [Hb], Mx ratio, Stutter percentage, etc...).

---

---

---

---

---

---

---

---

Some Semi-Continuous Examples

- LR mix (Haned and Gill)
- Balding (likeLTD - R program)
- FST (NYOCME, Mitchell *et al.*)
- Kelly *et al.* (University of Auckland, ESR)
- Lab Retriever (Lohmueller, Rudin and Inman)
- Armed Expert (NicheVision)
- Puch-Solis *et al.* (LiRa and LiRaHT)
- GenoProof Mixture (Qualitytype)

---

---

---

---

---

---

---

---

### Some Continuous Model Examples

- TrueAllele Casework (Cybergenetics)
- STRmix (ESR [NZ] and Australian collaboration)
- DNA-View Mixture Solution (Charles Brenner)
- DNAmixtures (Graversen 2013a,b) – open source, but requires HUGIN.

Weights may be determined by performing simulations of the data (Markov Chain Monte Carlo - MCMC).

---

---

---

---

---

---

---

---

### Goals of this Workshop

- To develop a greater understanding of the software systems presented.
- To foster discussions about training, validation, and scientific support of the software systems.
- To interact and ask questions of the software developers.

---

---

---

---

---

---

---

---

### Acknowledgments

National Institute of Justice and  
NIST OLES

Charlotte Word

Carol Bingham  
Lotte Downey



Points of view are mine and do not necessarily represent the official position or policies of the US Department of Justice or the National Institute of Standards and Technology.  
Certain commercial equipment, instruments, software and materials are identified in order to specify experimental procedures as completely as possible. In no case does such identification imply a recommendation or endorsement by the National Institute of Standards and Technology nor does it imply that any of the materials, instruments or equipment identified are necessarily the best available for the purpose.

Contact info:  
[mcoble@nist.gov](mailto:mcoble@nist.gov)  
**+1-301-975-4330**

---

---

---

---

---

---

---

---