

DTSA-II - QC

Nicholas W. M. Ritchie
Physicist, Microanalysis Research Group
National Institute of Standards and Technology
Gaithersburg, MD 20899-8371
nicholas.ritchie@nist.gov

Why?

- Is your instrument working correctly?
 - How can you demonstrate it?
- Was your instrument working correctly when a specific data set was collected?
 - How can you demonstrate it?
- Would you know if your instrument performance has degraded by 10%, 1%?
 - Are the standards you collected yesterday / a week ago / a year ago still good?


What could possibly go wrong?

- Bad energy scale calibration
 - Gain or zero offset out of calibration
- Noise at low energy or wacky Bremsstrahlung
 - Electrical, light, acoustic noise, ground loops
- Change in detector efficiency
 - Window contamination, detector position, obstruction, icing
- Deterioration in resolution
 - Electrical, light, acoustic noise, misconfiguration, ground loops
- Floating ground
 - Stage not grounded

Strategy

- Collect a spectrum
 - from a established material
 - at a consistent beam energy, probe current, live time
 - with a consistent working distance, geometry.
- Extract numeric metrics
 - Resolution, calibration, intensity etc
- Track the metrics
 - Plot them on QC charts

Example control charts



What does my QC look like?

- Sample: Pure Cu
- Beam energy: 20 keV
- Probe current: 1 nA
- Live time: 60 s
- Working distance: 17 mm
- Processor setting: Medium

I collect a new QC every day I collect a new set of quantitative microanalysis data.

950 × 301

Counts

15

10

5

16 18 20

Quality control alien

Previous: Welcome

Next: Select a measured spectrum

Select a QC project

EDS detector: Detector 2

QC Project: Pure copper at 20.0 keV on Detector 2

Material: Pure copper

Beam energy: 20.0 keV

Nominal working distance: 17.0 mm

Nominal probe current: 1.0 nA

☒ Use probe current normalization

Message: More...

Back Next Finish Cancel

Mass Fraction Atomic Fraction


Each detector can have one or more QC Project.

Using probe current normalization is a good idea!

The QC Project defines the sample and conditions under which the spectra are collected.

You define your own QC Project based on your own QC protocol.

You add data to a QC Project by collecting the appropriate spectrum and assigning it to the project.




Spectrum Report Command

Default Detector

MIRA3

Detector 2

Spectrum List



Previous: Select a QC project

Assign a measured spectrum

Next: Review the results

QC Spectrum

Spectrum

Probe current

Live time


Acquired

Message:

The date is read from the spectrum file and should be verified.

A QC is performed and the results are summarized and compared with the average, the earliest 10 measurements and the last 10 measurements.

When you select "Finish", the results are written to a database which can be used to generate comprehensive reports.



Spectrum Report Command

Default Detector

MIRA3

Detector 2

Spectrum List

- Pure copper[Wed Apr 23 080417 2014][2]
- Pure copper[Wed Apr 23 080417 2014][2]
- Fit[Pure copper[Wed Apr 23 080417 2014][2]]

Name	QC	Value	Average(All)	Avg(First 10)	Avg>Last 10)
Brem Counts		17616.607	17248.179±...	15162.162±...	17191.610±...
Channel width		9.9991±0.0...	9.9998±0.0...	9.9977±0.0...	9.9958±0.0...
Cu Ka		8783.620±1...	8248.026±4...	7501.397±1...	8521.800±1...
Cu Kβ		1228.574±5...	1144.741±5...	1022.609±1...	1186.614±...
Cu L-family		19161.918±...	18261.616±...	16329.023±...	18544.781±...
Dose		43.416	32.860±18....	58.579±3.462	43.204±2....
Duane-Hunt		19.8988	19.8649±0....	19.9247±0....	19.9065±0....
FWHM @ M...		132.476±0....	132.726±3....	131.836±0....	132.119±0....
Fano Factor		0.120	0.120±0.001	0.117±0.001	0.120
Noise		6.664±0.011	6.696±0.743	6.800±0.023	6.566±0.093
Total counts		2031524.000	1372116.53...	2344753.90...	1958776.10...
Zero offset		0.747±0.048	1.403±2.720	2.488±0.204	0.788±0.142

Message:

More...



Back

Next


Finish

Cancel

The QC spectrum and the fits are added to the Spectrum List.



A summary of the results is added to the daily Report.



QC Measurement Recorded

Spectrum	Pure copper[Wed Apr 23 080417 2014][2]
Index	357
Timestamp	2014-04-23 08:04:00.172
Project	Pure copper at 20.0 keV on Detector 2

Name	Value	First 10	Last 10	All
Brem Counts	17616.607	15162.162±379.378	17191.610±308.112	17248.179±8509.634
Channel width	9.9991±0.0002	9.9977±0.0115	9.9958±0.0061	9.9998±0.0193
Cu Ka	8783.620±14.224	7501.397±135.078	8521.800±144.388	8248.026±4017.949
Cu Kβ	1228.574±5.320	1022.609±18.033	1186.614±24.227	1144.741±558.908
Cu L-family	19161.918±21.008	16329.023±684.818	18544.781±296.238	18261.616±9043.992
Dose	43.416	58.579±3.462	43.204±24.191	32.860±18.005
Duane-Hunt	19.8988	19.9247±0.0395	19.9065±0.1277	19.8649±0.1140
FWHM @ Mn Ka	132.476±0.040	131.836±0.287	132.119±0.340	132.726±3.097
Fano Factor	0.120	0.117±0.001	0.120	0.120±0.001
Noise	6.664±0.011	6.800±0.023	6.566±0.093	6.696±0.743
Total counts	2031524.000	2344753.900±164191.258	1958076.100±1081892.462	1372116.534±792048.818
Zero offset	0.747±0.048	2.488±0.204	0.788±0.142	1.403±2.720

QUALITY CONTROL REPORT - APR 23, 2014


NIST DTSA-II Version Halley 2014-04-04

EPQ Version Halley 2014-04-04


Operator Tescan

Item	Value
Detector	Detector 2 - FWHM[Mn Ka]=128.0 eV - initial
Beam Energy	20.0 keV
Material	Pure copper


Channel width




Cu Ka




Cu Kβ



Cu L-family



You can also generate a more comprehensive QC report with control charts.



Select which data items you wish to summarize.

NIST DTSA-II

File Process Tools Report Help

Scaled to 60 nA-second

120,000 950 x 101

Pure copper[Wed Apr 23 080417 2014][2]
Fit[Pure copper[Wed Apr 23 080417 2014][2],Pure copper]+Bremsstrahlung[Pure copper
Pure copper[Wed Apr 23 080417 2014][2]-Fit[Pure copper[Wed Apr 23 080417 2014][2],F

Quality control alien

Previous: Select a QC project

Configure a report

Finish

Report items

- ☐ Brem Counts
- ☒ Channel width
- ☒ Cu Ka
- ☒ Cu K β
- ☒ Cu L-family
- ☐ Dose
- ☒ Duane-Hunt
- ☒ FWHM @ Mn Ka

Message: More...

Back Next Finish Cancel

Detector type Silicon Drift Detector
Detector window Moxtek AP 3.3 (manufacturer's table)
Elevation 35°
Energy Quadratic 0.00000000 eV/ch²
Energy offset 2.7 eV
Energy scale 10 eV/channel
Identifier Pure copper[Wed Apr 23 080417 2014][2]...
Instrument MIRA3

3 spectra selected

Mass Fraction Atomic Fraction

T

