


Standards vs References

Why do we sometimes need reference spectra in addition to standards?




Simple standards:

- Mg: Pure Mg
- Al: Pure Al
- Si: Pure Si
- Ca: CaF_2
- Fe: Pure Fe
- O: SiO_2

Like standards:

- Mg: K411 glass
- Al: Pure Al
- Si: K411 glass
- Ca: K411 glass
- Fe: K411 glass
- O: K411 glass



IUPAC	Standard	ZAF	Standard	ZAF
O All	SiO ₂	0.8037	K411	1.002
Mg All	Pure Mg	0.6695	K411	1.0325
Al All	Pure Al	0.6766	Al	0.6766
Si All	Pure Si	0.7569	K411	0.9558
Ca L-family	CaF ₂	0.43	K411	0.9873
Ca K-family	CaF ₂	0.9295	K411	0.9941
Fe L-family	Pure Fe	0.5981	K411	1.0031
Fe K α	Pure Fe	0.8195	K411	0.9948
Fe K β	Pure Fe	0.8204	K411	0.9947




Significant
corrections.



Close to
unity!


What is the problem?



These peaks are too close to each other.

Solution: Provide additional spectra which provide clear unobstructed perspectives of the obstructed peaks.

Solution



Fit CaF2 as a Ca reference, MgO as an O reference and Fe as a Fe reference to the K411 spectrum and to the unknown.

$$k_{unk,std} = \frac{k_{unk,ref}}{k_{std,ref}}$$