


# Metrology for Tissue Engineering

## Test Patterns and




## Cell Function Indicators

Scott Kennedy, Amit Sehgal, Newell Washburn (MSEL)  
John Elliott, Kurt Langenbach, Anne Plant (CSTL)  
Kenneth Yamada (NIH)  
Peter Jones (Colorado)

# Tissue Engineering: Defining Regeneration

## THE NEW ERA OF REGENERATIVE MEDICINE

*Dozens of biotech companies and university labs are developing ways to replace or regenerate failed body parts. Here are a few of the projects:*



### BONE

Bone-growth factors or stem cells are inserted into a porous material cut to a specific shape, creating new jaws or limbs. A product that creates shinbones is in clinical trials.

**COMPANIES:** Creative Biomolecules, Orquest, Sulzer Orthopedics Biologics, Genetics Institute, Osiris Therapeutics, Regeneron.



### SKIN

Organogenesis' Apligraf, a human-skin equivalent, is the first engineered body part to win FDA approval, initially for leg

ulcers. Other skins are in the works for foot ulcers and burns.


**COMPANIES:** Organogenesis, Advanced Tissue Sciences, Integra LifeSciences, LifeCell, Ortec International.



### PANCREAS

Insulin-manufacturing cells are harvested from pigs, encapsulated in membranes, and injected into the abdomen. The method has been tested in animals and could be in human trials in two years.

**COMPANIES:** BioHybrid Technologies, Neocrin, Circe Biomedical



### HEART VALVES, ARTERIES, AND VEINS

A 10-year initiative to build a heart has just started. Genetically engineered proteins have been successfully used to regrow blood vessels.

**COMPANIES:** Organogenesis, Advanced Tissue Sciences, Genetech, LifeCell, Reprogenesis.

DATA: BUSINESS WEEK, DRUG & MARKET DEVELOPMENT REPORTS



### SALIVA GLANDS

Proteins called aquaporins that allow cells to secrete water are used to recreate saliva glands damaged by disease or radiation. Glands are also being engineered to secrete healing drugs. The technique has proven successful in mice.

**COMPANIES:** None yet.



### URINARY TRACT

Cartilage cells are taken from the patient, packed into a tiny matrix, and injected into the weakened ureter, where they bulk up the tissue walls to prevent urinary backup and incontinence. The method is in late-phase clinical trials.

**COMPANIES:** Reprogenesis, Integra LifeSciences.



### BLADDER

Doctors at Children's Hospital in Boston have grown bladders from skin cells and implanted them in sheep.

They are about to try the same process on a patient.

**COMPANIES:** Reprogenesis.



### CARTILAGE

A product is already on the market that regrows knee cartilage. A chest has been grown for a boy and a human

**COMPANIES:** Genzyme Tissue, Biomatrix, Integra LifeSciences, Advanced Tissue Sciences, ReGen Biologics, Osiris Therapeutics



### TEETH

Enamel matrix proteins are used to fill cavities. It works in dogs; human trials are a few years away.

**COMPANIES:** Biora, Atrix Laboratories, Creative BioMolecules.



### BREAST

In preclinical studies, several companies have been able to create a cosmetic nipple by inserting a ball of cartilage. Researchers are now trying to grow a whole cosmetic breast.

**COMPANIES:** Reprogenesis, Integra LifeSciences.



### LIVER

A spongy membrane is built up and then seeded with liver cells. Organs the size of a dime have been grown, but a full-size liver could take 10 years due to its complexity.

**COMPANIES:** Advanced Tissue Sciences, Human Organ Sciences, Organogenesis.



### SPINAL CORD NERVES

Scientists are investigating nerve-growth factors, injecting them at the site of damage to encourage regeneration or seeding them along biodegradable filaments and implanting them. Rats have been made to walk again.


**COMPANIES:** Acorda, Regeneron, CytoTherapeutics, Guilford Pharmaceuticals.

*Generating functional replacements for missing or damaged tissue*

# Engineering Tissue by Controlling Cellular Response


Cellular responses characterized by:

- Proliferation
- Differentiation
- Matrix Production
- Migration
- Programmed Cell Death



Tissue engineering tools:

- Scaffolds
- Extracellular Matrix Proteins
- Growth Factors
- Gene Delivery



*Achieving this goal requires understanding complex signals that drive cellular differentiation and organization*

# Tissue Engineering: Current Products

## Skin substitutes

Apligraf (Organogenesis/Novartis)

Dermagraft (Smith-Nephew)

Dermal Regeneration Template (Integra Life Sciences)

## Cartilage engineering

Carticel (Genzyme)

Synvisc (Genzyme)

## Bone engineering

Smartscrew (Bionx Implants)

CMB401 (Integra Life Science/American Home Products)

## Dental tissue engineering

Biomend (Integra Life Sciences)


Over 70 start-up companies total\*

\*[Lysaght MJ, Reyes J. Tiss Eng 7, 485-493 (2001)]



*Industry needs a measurement infrastructure*






# The FDA Response: Office of Cellular, Tissue, and Gene Therapies

From the OCTGT Functional Statement:


[The Office of Cellular, Tissue, and Gene Therapies] plans and conducts research related to the development, manufacture, testing, and activities of cellular, tissue, and gene therapy products, including those related to AIDS and those prepared by genetic engineering and synthetic procedures, in order to develop and maintain a scientific base for **establishing standards designed to ensure the continued safety, purity, potency and effectiveness of cellular, tissue, and gene therapy products.**

# Metrology for Tissue Engineering



*Couple materials science and cell biology expertise to prepare reference materials of varying surface properties, cells that signal their functional state, and maps of cell responses*

# Measurements at the Bio/Material Interface




# Correlating Materials and Biology: The Deliverables

## Materials

- Polymeric reference materials to replicate real processing variables and direct cell response
- Reference biomimetic materials to validate indicator cell response


## Indicator Cells

- Develop genetically engineered cells for use as reporters of biocompatibility
- Demonstrate use of indicator cells of different tissue types for rapid screening of biomaterials




## Test Patterns and Cell Function Indicators

- Correlate cell response to materials properties with improved measurement methodologies






# Gradient Libraries: Spanning Material Parameter Spaces



# Metrology Chain for Gradient Libraries




Using gradient libraries that model manufactured materials or control cellular response:


- Validate processing techniques.
- Test uncharacterized cells and measure response.
- Calibrate indicator-cell response.

*Gradient library methodology with defined cell culture conditions can be used to validate cells and materials as a part of a complete metrology chain for tissue engineering.*


# Biomaterial Gradient Libraries




Polymer Crystallinity  
(nm roughness)




Polymer Blend Morphology  
( $\mu\text{m}$  roughness)




Surface patterns  
( $\mu\text{m}$  features)



Collagen Organization  
(ECM model)



Hydrogel Modulus  
(cellular mechanics)




Signaling-Peptide Density  
(nm organization)



*The performance of tissue engineered medical products is determined by both materials and processing conditions. Gradient libraries are capable of modeling both these parameters.*



# Automated Microscopy





- Spatially-addressable microscopy to probe two-dimensional array
- Fluorescent staining used to selectively image different cellular components
- Image analysis software for quantification of cellular response







# Cell Proliferation on Crystallinity Gradients



- Crystallinity gradient poly(L-lactic acid) film produced using temperature gradient.
- Osteoblasts cultured for 1, 3, and 5 d.
- Proliferation measured with automatic fluorescence microscopy.

# Quantitative Proliferation Data from Crystallinity Gradient Library




*Gradient libraries provide the dependence of cellular response as a function of material parameters.*

# Tissue Engineering Scaffolds



## Scaffold Parameter Space

- Polymer mechanical properties
- Polymer surface properties
- Polymer degradation kinetics
- Local geometry
- Characteristic length scales
- Anisotropy



*Methodologies are needed for measuring cellular response to scaffold structure and material properties.*


# Optical Coherence Tomography




*OCT is capable of providing three-dimensional structural information at micrometer-resolution, even in highly scattering media.*



# Functional Imaging: Coupling Fluorescence and OCT



OCT



Fluorescence image of  
stained cells in scaffold

*Simultaneous maps of scaffold structure and tissue development.*


# Correlating Materials and Biology: The Deliverables

## Materials

- Polymeric reference materials to replicate real processing variables and direct cell response
- Reference biomimetic materials to validate indicator cell response


## Indicator Cells

- Develop genetically engineered cells for use as reporters of biocompatibility
- Demonstrate use of indicator cells of different tissue types for rapid screening of biomaterials



## Test Patterns and Cell Function Indicators


- Correlate cell response to materials properties with improved measurement methodologies



# Systems Analysis Approach to Cell Biology

## Electronics


### The Black Box



## Systems Biology

### The Indicator Cell


### Biomarkers



***Cell-based indicators signal the result of a response to microenvironment.***

# Indicator Cells for Tissue Engineering

- Indicator Cells will provide a metric for rapidly assessing:
  - Cell response to new materials
  - Effect of materials processing/handling
  - Shelf life of materials
  - In vitro/In vivo biological response: inflammation, cell growth, biocompatibility, etc.
  - Cell response in 3-D matrices




***We will assist the development of TE products by providing indicator cells that allow real time evaluation of whether cells will proliferate, differentiate or apoptose when exposed to a new biomaterial.***



# Metrology Chain for Indicator Cells

## The Metrology Chain



### Using biomarkers that report cell response:

- Verify indicator-cell response.
- Calibrate response on reproducible defined materials.
- Test unknown materials and validate response on known materials.

***We will establish a metrology chain: NIST-validated markers used with NIST-validated conditions allows reliable assessment of new biomaterials.***

# Designing Indicator Cells

Identify genes, cell type  
Procure DNA  
Engineer artificial gene




•Verify DNA integrity: DNA sequencing, restriction enzyme analysis

Transfect cells  
Select clonal populations

•Verify DNA uptake:select for specific antibiotic resistance

Establish conditions for gene activation

•Confirm gene activation: Northern and Western blot analysis, immunohistochemistry



## Consultations on biomarkers

- ✎ Mark Pittenger, Osiris
- ✎ Glen Prestwich, Echelon
- ✎ J. Mansbridge, Advanced Tissue Sciences
- ✎ Peter Davies, IME UPenn
- ✎ Peter Jones, UCHSC
- ✎ R. Valentini, Cell Based Delivery


## Transfection vectors

- ✎ BD Biosciences / Clontech


## Verify signaling pathway

- ✎ Peter Jones, UCHSC


# Calibrating Indicator Cells



Reference  
surfaces of  
extracellular  
matrix (ECM)  
proteins




Quantitative  
fluorescence  
microscopy



Calibrate cell response  
on known reproducible  
surface

***Quantitative measurement of cell indicator response to reference conditions allows validation of indicator cells in other laboratories.***

# Indicator Cells in Practice: Smooth muscle cell proliferation in response to injury




Muscularized Artery

Injury



Blocked Artery from Hypertension




*Courtesy of Peter Jones, UCHSC*

Need for SMC Proliferation Control in TE

- Arterial Implants
- Synthetic Blood Vessels

# Smooth muscle cell proliferation in response to injury: Model System





*Several criteria are evaluated and used to validate the appropriateness and reliability of an indicator gene.*




# Reference Biomimetic Collagen Surfaces based on Self-Assembly



Denatured Collagen Film  
Thickness-  $3.6 \pm 0.2$  nm



Native Collagen Film  
Thickness-  $23 \pm 1$  nm



## WHY THIN FILMS based on self-assembly ?



- Currently no standardization of ECM surfaces.
- Removes surface inhomogeneity as a variable.
- Quantitative and reproducible.
- Could be fabricated in other laboratories.

***Reference materials will permit better inter-laboratory comparisons of data.***


# Validating Indicator Cell Response: GFP-Tenascin Expression

Assume that cell fluorescence  $\propto$  GFP expression  $\propto$  TN expression

Indicator Cell Status: Off  
Non-Proliferative



Indicator Cell Status: ON  
Proliferative




Optimization: Isolation of brighter single cell clones  
Degenerate GFP  
Transient transfection

# Validating Indicator Cell Response: Quantitative Microscopy



Fluorescence  
Microscope


X-Y sample stage



Objective  
Lens

Command

Image Analysis





Data Output

Cell Area: 3200


*Automated quantitative microscopy allows unbiased collection of data from a statistically relevant number of cells.*




# Validating Indicator Cell Response: Morphology and Proliferation



The mean size of SMCs, and the distribution of sizes, is highly reproducible on reference surfaces.





Proliferation rate is controlled by  
Reference Surfaces



**Morphology and proliferation measurements validate indicator cell response**

# Statistics and Population Distributions

## Cell Area on Collagen Films



Are these Distributions Similar?

	K-S test	$\chi^2$ test	T-test	F-test
Red & Blue	>95%	No	No	>95%
Red & Green	>95%	No	No	No
Blue & Green	>95%	No	>95%	No

K-S = Kolmogorov-Smirnov

**Identifying differences between population distributions may require rigorous statistical analysis**




# Correlating Materials and Biology: The Deliverables

## Materials

- Polymeric reference materials to replicate real processing variables and direct cell response
- Reference extracellular matrix materials to validate indicator cell response


## Indicator Cells

- Develop genetically engineered cells for use as reporters of biocompatibility
- Demonstrate use of indicator cells of different tissue types for rapid screening of biomaterials



## Test Patterns and Cell Function Indicators

- Correlate cell response to materials properties with improved measurement methodologies



# Accomplishments and Next Steps

## Accomplishments

- Development of gradient libraries that impact cell proliferation.
- Construction and preliminary validation of GFP-Tenascin indicator
- Molecular biology expertise and facilities
- Peptide synthesizer set up to synthesize signaling molecules
- Microscopy methods for quantitative analysis of cellular response

## Next steps

- Examine indicator cell response to polymer gradient materials
- Optimize TN promoter-GFP indicator cells
- Develop new promoter systems (for inflammation, oxidative damage, matrix production)
- Incorporate peptides into polymers and hydrogels to control cellular response
- Live cell quantitative microscopy

# Accomplishments and Next Steps

## Accomplishments

- Various cell lines established (bone, skeletal & smooth muscle, skin, kidney)
- Statistical methods developed for analyzing cell populations
- Established three-dimensional imaging methodologies.
- Construction of cell biology laboratories
- Reference surfaces of collagen
- Nine manuscripts in press, in preparation or submitted; 15 presentations made, 10 invited

## Next steps

- Develop transfection protocols for other cell types
- Validate utility of parametric and non-parametric fitting
- Characterize engineered bone and cartilage
- Investigation of broader range of cellular responses
- Develop new ECM reference surfaces

# Building the Infrastructure

## *Polymers Division infrastructure*

- Renovated 2000 ft<sup>2</sup> for biomaterial-library preparation, cell biology and imaging laboratories.
- Expand efforts through hiring staff and post-doctoral associates.

## *Biotechnology Division infrastructure*

- Capitalize on existing thin film fabrication and analysis capabilities.
- Expand capabilities and expertise in molecular and cell biology by hiring Kurt Langenbach.

# Strategic Partnerships

## *Establish strategic partnerships outside of NIST*

Consultation on Gradient Libraries: ●NIH (Kenneth Yamada), ●UPenn (Paul Jamney), ●Integra Life Sciences (Fred Cahn), ●NIH (Rocky Tuan)

Consultation on Biomarkers: ●Osiris (Mark Pittenger), ●Echelon (Glen Prestwich), ●Advanced Tissue Sciences (J. Mansbridge), ●IME, UPenn (Peter Davies), ●UCHSC (Peter Jones)

## Other collaborations

- Robert Levi, CHOP (chemically modified stents); ● Joachim Kohn, Rutgers (Preparation of new gradient libraries); ● Richard Spencer, NIH (Characterization of engineered cartilage); ● Clontech (Specific, high efficiency vectors); ● ReGen Corp (Collagen structure and function); ● Richard Gross, Brooklyn Polytechnic (High-throughput investigation of protein adsorption); ● Becton-Dickinson (Extracellular matrix substrates for comparison testing); ● Leica (Evaluation of serum protein interaction with materials); ASTM Workshop on Metrology for Cell Signaling ● Development of reference scaffolds with ASTM involvement

*By working with strategic partners we will provide measurement tools that industry can use to qualify materials, cells, and manufacturing practices*